

When Hearing Aids are Not Enough

Overall satisfaction with hearing instruments is high, nevertheless 95% of hearing instrument wearers want better understanding in noise. An FM system is the solution!

An FM system is a radio transmission system that consists of 2 parts:

1) **Transmitter:** has a microphone which picks up the desired sound and sends it on by radio waves to the receiver

2) **Receiver:** receives the radio signal from the transmitter and sends it to the hearing instruments

An FM system means that hearing over distance, in noise and in reverberant environments is made easy!



Phonak SmartLink Transmitter



Widex Scola Transmitter



Phonak MyLink Receiver



Manningham Hearing Clinic

...making hearing possible

Ian C Skipworth MACAud
m: 0418 377 875
e: hear@acousticaid.com.au
www.acousticaid.com.au

t: (03) 9857 3993
f: (03) 9857 4994
3/385 Belmore Road,
Balwyn VIC 3103