

# Open Fit Hearing Instruments


Most hearing aid manufacturers are now offering “Open Fit” hearing aids. What are the benefits of Open Fit aids?

**Open mould/dome:** Having a very open mould/dome in your ear canal makes a big difference. You will not feel like you are wearing hearing aids (except for the great sound quality)! The Occlusion effect (the blocked up sensation when something is in your ear) is not an issue with Open Fit hearing aids.

**Mini size:** Open Fit aids come in various sizes. Most of them are very small. These fit discreetly behind your ear and are very comfortable to wear. Some of the styles are illustrated below.

**Slim tubing system:** This slim tube system is practically invisible and is very lightweight to wear, especially when compared to the traditional BTE tubing.

Open Fit hearing aids are very easy to keep clean and free of wax with the cleaning tools provided at the time of fitting.



**Manningham Hearing Clinic**

*...making hearing possible*

**Ian C Skipworth MACAud**  
m: 0418 377 875  
e: [hear@acousticaid.com.au](mailto:hear@acousticaid.com.au)  
[www.acousticaid.com.au](http://www.acousticaid.com.au)

t: (03) 9857 3993  
f: (03) 9857 4994  
3/385 Belmore Road,  
Balwyn VIC 3103